

***SHRIMP CHECK***

detection kits .

虾类疾病检测试剂盒系列



**Proven shrimp diseases detection kits  
are now within reach**

全方位虾类疾病检测试剂盒系列

Rapid Test Kit . Extraction Kit . PCR Kit . qPCR Kit

快速检测试剂盒 . PCR试剂盒 . 实时荧光定量PCR试剂盒

# Early detection via ShrimpCheck may PROTECT your shrimp against HIGHLY INFECTIOUS DISEASES

ShrimpCheck 虾类疾病检测试剂盒系列，有助于虾苗选购和养殖管理，由源头避免虾苗感染

## ShrimpCheck Shrimp Diseases Detection Kits ShrimpCheck 虾类疾病检测试剂盒系列



**Rapid Test Kit**  
快速检测试剂盒



**PCR Kit**  
PCR 试剂盒



**Real-Time PCR Kit**  
实时荧光定量 PCR 试剂盒

<b>Types of Kits</b> 种类	WSSV, IMNV, YHV	WSSV, MBV, IMNV, YHV, EHP, AHPND/EMS, HPV, IHHNV, TSV, CMNV	WSSV, MBV, IMNV, YHV, EHP, AHPND/EMS, HPV, IHHNV, TSV, CMNV, SHIV/DIV1
<b>Assay Time</b> 检测时间	Approx. 15 minutes 约15分钟	Approx. 3-4 hours 约3-4小时	Approx. 2-3 hours 约2-3小时
<b>Skilled Operator</b> 专业操作人员	Not required 无需	Required 必需	Required 必需
<b>Features</b> 特点	On-site and rapid screening 现场快速检测	Accurate positive control to confirm findings 附上正控制组以确认结果	Accurate positive control to confirm findings 附上正控制组以确认结果
	Simple to use 操作简单	Internal/external control is provided, eliminates false negative results 附上内/外控制组，摒除假阴性结果	Internal/external control is provided, eliminates false negative results 附上内/外控制组，摒除假阴性结果
	No expensive tools required & cost effective 无需昂贵器材与耗材，经济实惠	High specificity and accuracy 高敏感度与准确性	High specificity and accuracy 高敏感度与准确性
			Easy interpretation of the results without gel electrophoresis 无需凝胶电泳步骤，轻易分析结果

**Speedy Assay Sdn. Bhd.** 200201018654 (586317-W)

### Headquarters

Revongen Corporation Center  
Level 17 • Top Glove Tower  
No. 16 • Persiaran Setia Dagang • Setia Alam • Seksyen U13  
40170 Shah Alam • Selangor Darul Ehsan • Malaysia.

t • +6 03 3359 1166  
f • +6 03 3358 0303  
e • info@speedyassay.com  
w • www.speedyassay.com



Follow Us

Speedy Assay